


Injamam - Boronic acid- Test tube 1- C nmr
C13CPD256 CDCl3 /opt/topspin Adv_Exp 41



NAME X_20160114
EXPNO 541
PROCNO 1
Date_ 20160115
Time 9.53
INSTRUM spect
PROBHD 5 mm BBO BB-1H
PULPROG zgpg30
TD 65536
SOLVENT CDCl3
NS 256
DS 4
SWH 23980.814 Hz
FIDRES 0.365918 Hz
AQ 1.3664756 sec
RG 18390.4
DW 20.850 usec
DE 6.50 usec
TE 303.0 K
D1 2.00000000 sec
D11 0.03000000 sec
TD0 1

===== CHANNEL f1 =====
NUC1 13C
P1 8.90 usec
PL1 6.00 dB
SFO1 100.6479773 MHz

===== CHANNEL f2 =====
CPDPRG2 waltz16
NUC2 1H
PCPD2 80.00 usec
PL2 3.00 dB
PL12 21.00 dB
PL13 21.00 dB
SFO2 400.2316009 MHz
SI 32768
SF 100.6379140 MHz
WDW EM
SSB 0
LB 1.00 Hz
GB 0
PC 1.40