

B.DRAIN AND C>BEVAS PURE 8
Check purity before reaction start
PROTON CDCl3 /opt/topspin Adv_Exp 7



NAME X_20160224
EXPNO 10
PROCNO 1
Date_ 20160224
Time 9.48
INSTRUM spect
PROBHD 5 mm BBO BB-1H
PULPROG zg30
TD 65536
SOLVENT CDCl3
NS 16
DS 2
SWH 8278.146 Hz
FIDRES 0.126314 Hz
AQ 3.9584243 sec
RG 181
DW 60.400 usec
DE 6.50 usec
TE 303.0 K
D1 1.00000000 sec
TD0 1

===== CHANNEL f1 =====
NUC1 1H
P1 10.60 usec
PL1 3.00 dB
SFO1 400.2324716 MHz
SI 65536
SF 400.2300370 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00

